



VetLogic End of Life



AVIMARK®


End of Life for Vetlogic

VetLogic was a method to transmit and receive digital information. This method is being replaced with a new solution and VetLogic, via the Avimark Sync Service, must be disabled.

Avimark Sync Service

To disable VetLogic for Avimark, it is necessary to change the Avimark Sync Service from True to False.


1. In Avimark, click  to open Options Maintenance.
2. In the Search bar, type **sync** to display Avimark Sync Service.
3. Select **Active**.



4. Click **Change**.
5. Change True to **False**.
6. If a warning displays to confirm disabling the Sync Service, click **Yes**.
7. Click **Done** to close Options Maintenance.
8. Close all instances of Avimark on all computers until the remaining steps are complete.


VETLOGIC END OF LIFE

9. On the **server**, locate the shortcut to Avimark Guardian then **right-click** > **Run as Administrator**.



10. The Avimark Server Guardian window opens.

- a. Click **Stop the Service**.
- b. Click **Run**.
- c. Click **Start the Service**.
- d. Click **Run**.



11. Log back into Avimark.