



VetWeb Online Access Instructions



855-478-7920 | 304 Ohio St. | Oshkosh, WI 54902
Revised July 31, 2019

Table of Contents

Additional Instructions for Accessing VetWeb Online	3
Mac Users	3
Windows Users	4

Additional Instructions for Accessing VetWeb Online

The instructions below will help guide new clients through the steps to access VetWeb Online and their Covetrus Practice Management Software (PiMS). Instructions are available for both Mac Users and ["Windows Users"](#).

Mac Users

1. Download and install the **Microsoft Remote Desktop** app.
2. Open the Microsoft Remote Desktop app.
3. Click on **Preferences**.
4. Click on **Gateway**.
5. Click on the "plus" sign to add a Gateway.
6. Enter **VetWeb Online** as the **Gateway name** and **vetwebonline.com** as the **Server**.
7. Close Preferences.
8. Add the PiMS information:
 - a. Click **New**.
 - b. **Connection name**: type **VetWeb Online**.
 - c. **PC name**, enter one of the following:
 - For AVImark: **aviproduct.vetwebonline.com**
 - For ImproMed: **impprod1810.vetwebonline.com**
 - For ImproMed Equine: **impeqprod1810.vetwebonline.com**
9. For the **Gateway**, select **VetWeb Online**.
10. For **Credentials**, enter the user name and password assigned to you by your Conversions Specialist.
11. For **Resolution**, **Colors**, and **Full screen mode**, leave the default settings.
12. Check the box labeled **Scale content**. The options, *Start session in full screen* and *Use all monitors* should already be checked.
13. Close the setup screen.
14. From **My Desktop**, double-click on **VetWeb Online** to run.

Windows Users

To open VetWeb Online on a Windows-based PC:

1. Click **Start | Run** or **Start | Search**.
2. Type **mstsc** then press **Enter**.
3. Click **Show Options**. *This expands the dialog box to view more General settings and more connection tabs.*
4. In the **Computer** field, type one of the following to match your PiMS:
 - o For AVImark: **aviproduct.vetwebonline.com**
 - o For ImproMed: **impprod1810.vetwebonline.com**
 - o For ImproMed Equine: **impeqprod1810.vetwebonline.com**
5. Enter the **User name** provided by your Conversions Specialist. *If you do not have a User name, contact a member of the Systems Implementations team.*
6. Click on the **Advanced** tab.
7. Click on the **Settings...** button.
8. Choose the option, **Use these RD Gateway server settings**.
9. Enter **vetwebonline.com** for the **Server name**.
10. Select the option, **Allow me to select later** as the **Logon method**.
11. Check the option, **Bypass RD Gateway server for local addresses**.
12. Check the option, **Use my RD Gateway credentials for the remote computer**.
13. Click **OK** to close Settings.
14. Click **Connect** to open VetWeb Online.

If you continue to have difficulty connecting to VetWeb Online, please contact a Systems Implementations team member at **(855) 478-7920**.